lospital ID:	BLOOD STORAGE UNIT FAILURE
	TEMPERATURE DOCUMENTATION

FAX DAILY TO MEDIC.

If blood is not stored in an approved blood storage unit, the temperature must be monitored and recorded frequently. Determine frequency based on the reliability of the alternate storage unit, but MEDIC recommends no less frequently than every 2 hours. If the alternate storage has been validated to maintain proper temperature for a longer period, a maximum interval of recording temperature every 4 hours could be used.

Date:			Date:				
Time	Temp	Tech		Time	Temp	Tech	
11110	Temp	Teen		111110	Temp	10011	
Comments:							
Comments.							
Supervisory Review By:							
			Date:				
MEDIC Review By:			Date:				
MEDIC Regional Blood Center							