

Attach Recorder Chart Here

## DAILY TEMPERATURE DOCUMENTATION FOR RED BLOOD CELL STORAGE

Hospital ID: \_\_\_\_\_

Blood Storage Unit: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

\_\_\_\_\_

Recorder Chart Changed: Day: \_\_\_\_\_ Date: \_\_\_\_\_ Time: \_\_\_\_\_ By: \_\_\_\_\_

Day	Date	Time	Recorder Temperature °C	Manual/Digital Thermometer °C		A/U	Tech	Check Inventory Level
				Upper	Lower			

**NOTE:** Acceptable storage temperature for Red Blood Cells is 1 to 6 °C.  
Recorder chart and upper and lower thermometers must read within 1 °C.

**KEY:** A = Acceptable      U = Unacceptable

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Transfusion Service Supervisory Review By: \_\_\_\_\_ Date: \_\_\_\_\_