

REFERENCE LABORATORY

I. REFERENCE LABORATORY HOURS OF OPERATION:

Regular Hours - Monday through Friday 8:00 AM to 4:00 PM.

On-Call Staff - Evenings, Nights, Weekend and Holidays.

NOTE: *All Reference Lab procedures are subject to after-hours charge if requested to be performed outside regular hours.*

If there is not an urgent need for blood, the work up will be done by the Reference Lab during the next normal working day (Monday – Friday 8:00 AM to 4:00 PM).

”Urgent” is defined as a clinical situation where there is active bleeding or unstable anemia with Hgb/Hct levels at or approaching symptomatic or life threatening levels. Blood to be available for elective surgery or an outpatient is not defined as urgent.

II. HOW TO REQUEST REFERENCE LAB ASSISTANCE:

A. CALL REFERENCE LAB:

Call (865) 521-2696 and request to speak with a Reference Lab Technologist before sending any sample. Reference Lab will ask for:

1. Patient’s name:

Has MEDIC seen this patient before?

When?

2. Serological Problem:

What problem have you encountered?

Do you have an ABO discrepancy, positive antibody screen, incompatible crossmatch, positive DAT, transfusion reaction, etc.?

3. Crossmatched Units:

Have you crossmatched any units?

If yes, how many?

How many were compatible?

How many were incompatible?

4. Blood Needs:

Does the patient need blood?

How many units? By when?

Is patient bleeding?

REFERENCE LABORATORY

What is the patient's hemoglobin?
Is patient scheduled for surgery?

5. Patient's Diagnosis and History:

Has patient ever been pregnant or transfused?
When was the last transfusion?
(Be sure to check blood bank records.)
Has the patient received Darzalex?

6. Previous Antibodies Identified:

Does the patient have a history of an antibody identified by any other consultation blood bank or service?

B. REQUEST FOR SEROLOGICAL PROBLEM RESOLUTION:

Complete Request for Serological Problem Resolution:

NOTE: *Failure to provide all requested information may result in erroneous testing.*

C. SAMPLE REQUIREMENTS:

Send adequate fresh sample (less than 24 hours old) when submitting a specimen for serological problem resolution. Additional tubes will be requested if insufficient sample is sent or if the serological problem is such that more serum or cells are necessary. All tubes should be properly labeled with patient's name, a second unique identifier, and date the sample was drawn.

MINIMUM: *One full 7 mL EDTA Tube
One full 10 mL Red Top Tube (no serum separators)*

NOTE: *Failure to send the minimum sample volume may result in the work-up being delayed or not completed.*

D. TRANSPORTATION OF SAMPLE TO MEDIC:

MEDIC can arrange pick-up of specimens for Reference Lab workups upon request, or specimens may be delivered to MEDIC's Hospital Services loading dock area. Clear communication is key!

REFERENCE LABORATORY

III. RECEIPT OF SAMPLE BY REFERENCE LAB:

A. ROUTINE SEROLOGICAL PROBLEM WORKUP:

The routine serological problem workup will include:

1. ABO and Rh
2. Selected Phenotype (C, E, c, K and e if needed)
3. DAT
4. Antibody Screen and Identification (if applicable)
5. Patient Testing (if requested) with antigen negative units

NOTE: *Units tested with the patient's sample must be crossmatched at Transfusion Service prior to issue.*

B. COMPLICATED SEROLOGICAL PROBLEM RESOLUTION:

Additional tests performed to resolve a more complicated serological problem will depend upon the individual case and the patient's history.

C. REPORTS:

1. Telephone Report:
A telephone report will be made when serologic tests are completed.
2. Written Report:
A written preliminary report will be mailed/faxed or sent with the requested units.
3. Reference Lab Report Summary:
A written final report will be mailed/faxed upon completion and review of all results.

IV. CHANGE IN PATIENT STATUS:

Communication is essential. Please contact MEDIC and the receiving hospital if the patient is transferred to another hospital or the physician changes or cancels the request.

V. REQUESTS FOR ANTIGEN NEGATIVE UNITS:

All after hour requests for antigen negative units should be made through BloodHub.

NOTE: *ABO compatible units will be issued after regular Reference Lab hours.*

REFERENCE LABORATORY

VI. REQUESTS FOR PLATELETS:

PLATELET CROSSMATCHING:

Platelet crossmatching is available Monday - Friday from 8:00 AM to 4:00 PM. Call the Reference Lab at (865) 521-2696 before sending sample. We require one freshly drawn patient sample (plasma or serum-no gel tubes) and request sample delivery to MEDIC by 1:00 PM for same day testing. Deliver sample and completed Request for Platelet Crossmatch to the Hospital Services Department.

VII. REQUESTS FOR HEMOGLOBIN S NEGATIVE RBC PRODUCTS:

Hemoglobin S screening is available Monday – Friday from 8:00 AM to 4:00 PM. Call the Reference Lab at (865) 521-2696 to request.

NOTE: *ABO compatible units will be issued after regular Reference Lab hours.*